İmalatçı:

OLİGA Mak. Elekt. İth. İhr. San. ve Tic. Ltd. Şti.
İvedik OSB Mah.1453.Cad. 5/D Yenimahalle/Ankara/TÜRKİYE
Tel: +90 312 220 15 78
Faks: +90 312 220 15 76
www.oliga.com.tr
info@oliga.com.tr

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Information for Use
Please read prior to use

INTRAMEDULLARY NAILING SYSTEMS

(Refer to product catalog for product technical details, names and reference number.)

RESTRICTIONS FOR REUSE

Intramedullary Nailing Systems are designed for SINGLE USE. Do not reuse the devices.

PRODUCT CONTENT AND CONFORMITY

The implantable parts of Intramedullary Nail Systems are manufactured out of titanium. (Ti-6Al-4V). One of the nails' internal components, the Compression Bolt Screw, is manufactured from cobalt-chromium-molybdenum (Co-Cr-Mo). The instruments are manufactured from stainless steel, aluminum or carbon fiber reinforced PEEK.

Products do not contain human blood derivatives, animal tissues, drugs and phthalates.

Literature studies have shown that titanium and cobalt-chromium-molybdenum products are MR-compatible. Intramedullary Nailing Systems were not subjected to any specific test to evaluate the safety and compatibility with the MR environment.

Therefore, OLİGA Mak. Elekt. İth. İhr. San. ve Tic. Ltd. Şti. cannot make any recommendations regarding the safety and image accuracy in MR procedures. Nails have not been evaluated for safety and compatibility in the MR environment.

DESCRIPTION

Intramedullary Nailing Systems are designed for fractures of the femur and tibia. Depending on the preferred nail model, the system consists of proximal and distal end cap, locking screws, bolt screws and pins. These implants are not designed to replace the normal healthy bone. The loads produced from patient's weight bearing and the activity levels will influence on the implant's life. The patient should be informed that stress on the implant can be produced even if there is no load on the limb. Without robust bone union, the limb's weight, muscle strength associated with a moving limb, or repeated small stresses can cause the implant to fail. For this reason, the patient must comply fully with the post-operative instructions given by the surgeon.

Intramedullary Nailing Systems should be used only by surgeons having full knowledge of orthopedic fixation procedures. Surgeons must be familiar with the use of implants and instruments before surgery, including removal procedures. In order to obtain the Surgical Technique file containing detailed surgical information, please contact OLİGA Mak. Elekt. İth. İhr. San. ve Tic. Ltd. Şti. or your local supplier.

INDICATIONS

Intramedullary Nailing Systems are intended to be implanted in the medullary canal of long bones, femur and tibia) for the treatment of:

- Trauma fractures
- Refractures
- Nonunions/Malunions
- Malalignments,
- Pathological fractures and
- Bone deformities

CONTRAINDICATIONS

The contraindications below apply to all product categories.

The following conditions may increase the risk of implant failure. These conditions cannot be presented in full, as they may vary depending on the patient and adverse events that may develop during or after the operation. Doctors should use their clinical experience while deciding which implant should be appropriate for a particular patient.

1. Local active infection.
2. General medical conditions including defective blood supply, pulmonary insufficiency or latent infection.
3. Patients who do not want or cannot follow the instructions of the surgeon after the operation.
4. Obese patients.
5. Patients with documented metal allergy/intolerance or suspicion.
6. Patients with obliterated medullary canals.
7. Comminuted fractures.
8. Patients with insufficient bone quality to support implants.
9. Patients with cognitive and/or physical disabilities that may lead to unacceptable risks of failure of fixation failure.

GENERAL WARNINGS

1. Implants and their parts are designed for SINGLE USE and should NEVER be reused.
2. It is known that the use of cigarettes, chronic use of steroids and other anti-inflammatory drugs may affect bone healing and have a negative effect on bone regeneration.
3. Care should be taken to use and store the products. Cuts, bends and scratches may significantly decrease the functional strength of the components
4. None of the intramedullary nails and their parts have been approved for use in any area of the body except for the ones to which they are intended.
5. When using Intramedullary Nailing Systems, do not use products and instruments from other manufacturers as they may be manufactured from different materials and cause adverse reactions.

6. All products should be carefully examined before use. Product integrity, sterility and performance are only provided if the product packaging is not damaged. If the packaging is damaged or a part of the product is damaged, DO NOT USE the product.
7. Instruments should be kept away from abrasive media such as saline water while storing. Before the operation, verify the intactness and undamaged state of the implants by conducting a control.
8. Selecting the appropriate implant model and size is critical. Larger and longer nails may be needed for adequate stabilization of severely fractured or unstable regions. Use of an implant at inappropriate size or improper positioning of the implant can cause loosening, bending, cracking and fracture of the implant within the bone, or a combination of those of those.
9. The length of the locking screws shall be measured by the method specified in the Surgical Technique file.
10. When the locking screws of Intramedullary Nailing Systems are fitted, care must be taken that the screws do not protrude beyond the joint surface or excessive protrusion into soft tissue.
11. The nails and lag screw pins can only be deployed as described in the nail's surgical technique using the appropriate tools. The use of a different tool or a powered device can damage the system and prevent good fixation or implant removal.
12. It is important that the pins or locking screws are not placed close to the fractured region. This can adversely affect fixation.
13. The canal should be enlarged depending on the selected nail diameter, type of fracture, bone quality and diameter.
14. If the implant cannot be inserted during surgery, the surgeon should never continue to hammer or increase force. If necessary, the implant should be removed and the medullary canal should be enlarged or a smaller implant should be considered.
15. The length of the end cap must be carefully selected to avoid damage to the joint.
16. Radiologic evidence of callus formation is needed in order for the patient to be allowed full weight bearing full. Partial weight bearing is recommended when there is good contact between two intact parts of the bone that allows for load sharing. Hip and knee motility should always be encouraged within pain limits. However, good clinical outcomes are achieved with the earliest possible mobility and tolerable weight according to the patient's condition.
17. The patient's healing process should be carefully monitored. If fracture healing is slow, fracture healing may be promoted by other measures, such as implant dynamization, grafting, or implant replacement.
18. Any adverse or unexpected event should be recorded into surgeon reports.
19. The surgeon should decide whether or not to remove an implant. In general, intramedullary nails are removed about one year after bone healing, if necessary. After removal of the implant, adequate post-operative controls should be performed to prevent recurrence. The patient should be warned that a second operation, usually of a smaller duration, is required to remove the implant. Removal of the implant is not recommended unless there is pain in elderly or weak patients.
20. For at least 6 months following the operation, the patient should be x-rayed periodically to obtain information about the progression of the healing process.

RISKS OF REUSING A "SINGLE USE" PRODUCT

Intramedullary Nailing Systems are composed by reusable instruments and disposable implants. The disposable single use implants are identified by the symbol  on the product label.

Implant: A product intended to be completely/partially placed in the human body by surgical intervention and intended to remain in place for at least 30 days after the procedure can be considered an implant. Once removed from the patient, the implant should be destroyed. Reuse of the implant poses a risk of contamination for users and patients. Repeated use of the implant creates a risk of contamination for users and patients. Reuse of the implant cannot guarantee the mechanical and functional performance of a product that has never been used, jeopardizes the efficacy of the products and brings health risks for patients.

POSSIBLE SIDE EFFECTS

8. If incorrect loads are applied, they can cause loosening, twisting or breaking of implanted parts leading to malunions or nonunions.
9. Loss of anatomic position due to malunion.
10. Scar formation possibly causing pain and/or neurological compromises around nerves.
11. Intrinsic risks associated with anesthesia and surgery (hemorrhage, hematoma, seroma, embolism, edema, excessive bleeding, phlebitis, wound/wound infection, damage to blood vessels or nerves).
12. Nonunion or delayed union, which may lead to implant failure.
13. Metal sensitivity or allergic reaction to a foreign body.
14. Pain, discomfort or abnormal perception due to the presence of the implant.

IMPORTANT INFORMATION

A successful result cannot be achieved in every surgical case. Additional complications arising from improper use, medical reasons or product failure may develop. These additional complications may require surgical intervention to remove or replace the implant.

Knowledge of pre-operative procedures, including surgical procedures, appropriate product selection and placement, are important considerations in the successful use of implants.

The best performance of the pins is achieved by deploying to the cortical bone on the bone shaft. The pins should not be deployed in areas where the cortex is expanding (metaphysis). In these cases, a shorter implant should be chosen.

Adequate patient selection and the ability of the patient to adhere to the physician's instructions and post-operative instructions greatly influence the outcome. It is important to choose the most appropriate treatment according to the patient's physical/mental activity needs and limitations. If a patient is contraindicated or susceptible to any contraindication, DO NOT USE Intramedullary Nailing Systems.

Post-operative care is extremely important. If the patient does not follow the post-operative instructions, the implant may fail. The patient should be warned that revision surgery will be necessary to remove the implant in those cases.

STERILIZATION

All implantable parts of Oliga Intramedullary Nailing Systems are available for sale as STERILE. All instrument used with the products are provided as NON-STERILE. Non-sterile products require cleaning and sterilization prior to use. Refer to the products label to determine the sterility of the products.

PACKAGING AND LABELING

The parts should be only accepted by the hospital or by the surgeon when delivered in undamaged factory packaging and label. Labels in the product are for the customers and for the hospital.

PRODUCT DETAILS AND EXPLANATIONS

Intramedullary Nailing Systems are not designed for any purpose except for what is described below.

1 PROXIMAL FEMORAL NAIL SYSTEM

Nails with different lengths and lag targeting angles are available.

Indications:

- Intertrochanteric fractures
- Pertrochanteric fractures
- High subtrochanteric (without shaft extension)
- Patological fractures
- Malunion/Nonunion
- Combination of those fractures
- Long nails are indicated for low subtrochanteric fractures

Contraindications:

- Distal third fractures
- Femoral neck fractures.

2 FEMORAL NAIL SYSTEM

Nails with different diameters and lengths are available.

Indications:

- Femoral shaft fractures
- Ipsilateral hip/shaft fractures
- Ipsilateral femur/tibia fractures (floating knee)
- Supracondylar fractures including those with intraarticular extension
- Fractures proximal to a knee implant
- Fractures distal to a hip implant
- Osteoporotic fractures
- Pathologic fractures
- Malunions/nonunions

3 TIBIAL NAIL SYSTEM

Nails with different diameters and lengths are available.

Indications:

- Proximal tibial fractures
- Open and closed fractures of the tibial shaft
- Segmental fractures
- Pathological fractures
- Malunions/nonunions

Contraindications:

- Tibial Nails are not indicated for fractures in the distal tibia.

The usage symbols for the product are given on the previous page, Article "Symbols and Meanings".